



Guide to Cementation Selection





Variolink® II


ESTHETIC RESIN CEMENT

Features

- Proven long term clinical success
- Superior strength with outstanding esthetics
- (6) Light cure base shades and (2) catalyst shades
- (3) Viscosities
- Highest Radiopacity in its category



Variolink® II's excellent radiopacity permits the clear distinction of the restoration, cement, and caries on x-rays.




Appeal™

ESTHETIC RESIN CEMENT

Features

- Unique chemistry provides fast and easy clean up
- Highly translucent "Value" shades enhance the esthetics of the final restoration
- (7) Light Cure, (3) Dual cure values
- Microfill chemistry provides low wear rates and high polishability
- Amine free for outstanding color stability




Multilink®

UNIVERSAL ADHESIVE CEMENT

Features

- High bond strengths provide the ultimate security for all types of restorations.
- Universal applications make Multilink an ideal cement for indirect restorations (inlays, onlays, crowns, bridges, and posts) made of:
 - Metal and metal ceramics (PFM)
 - Metal-free (all-ceramics)
 - Strengthened core ceramics (Aluminum oxide, zirconia)
 - Composites
- Multilink applies quickly and cleans up easily.
- Self-etching and self-curing
- Three shades with distinct translucencies to meet all clinical indications.



Vivaglass® CEM

GLASS IONOMER CEMENT

The Easiest Method of Conventional Cementation

Features


- Convenient injection capsule
- Predosed mixing ratio
- Flexible applicator tip
- Low film thickness
- High radiopacity
- Continuous fluoride release
- Ideal capsule volume - enough for 1-2 crowns



Selecting the Right Resin Cement

Recent developments in restorative materials and cements continue to support practitioners' ability to deliver aesthetic dental care. The proliferation of these therapeutic alternatives can be challenging, particularly when the clinician is

treatment planning and selecting the proper techniques, ceramic systems, and cements. The following guide highlights the necessary steps and benefits associated with three cements that are currently available at Ivoclar Vivadent.



When the clinician has three quality cements from which to choose, as shown herein, the ideal selection for a given indication often amounts simply to a matter of personal preference. Multilink, Appeal, and Variolink II can all be used chairside for optimal aesthetic results.

(Photographs and dentistry courtesy of Chris Ramsey, DDS.)
 *Cast gold crown was also seated with Multilink.

Appeal™

ESTHETIC RESIN CEMENT

Appeal Light Cure Anterior Assortment – 574341

Appeal Light Cure Base Values - 1x1.5 ml

- 574340 High Value +3 (Bleach)
- 574348 High Value +2 (B)
- 574347 High Value +1 (Transparent)
- 574346 Medium Value 0 (Clear)
- 574344 Low Value -1 (A2)
- 574345 Low Value -2 (A3)
- 574342 Low Value -3 (A4)

Appeal Base Try-In Values - 1x1.5 ml

- 574332 High Value +3
- 574339 High Value +2
- 574338 High Value +1
- 574337 Medium Value 0
- 574334 Low Value -1
- 574336 Low Value -2
- 574333 Low Value -3

- 1 Monobond® S
- 1 Liquid Strip®
- 1 Total Etch®
- 1 Excite® Bottle
- 1 Excite® Applicators
- 1 Brush holder
- Various canulas and brushes
- 1 Vivapad
- 1 Mixing pad
- 1 Appeal Light Cure Shade Guide
- Instructions for use and flowchart



Appeal Dual Cure Posterior Assortment - 574343

Appeal Dual Cure Values – 1x3 ml

- 574330 High Value +1 (Transparent)
- 574331 Medium Value 0 (A1)
- 574329 Low Value -1 (A4)

- 1 Monobond® S
- 1 Liquid Strip®
- 1 Total Etch®
- 1 Excite® Bottle
- 1 Excite® Applicators
- 1 Brush holder
- Various canulas and brushes
- 1 Vivapad
- 1 Mixing pad
- 1 Appeal Dual Cure Shade Guide
- Instructions for use and flowchart



Variolink® II

ESTHETIC RESIN CEMENT

Esthetic Cementation System - 558946

Base – 1x2.5 gm

- 558949 Bleach XL
- 558950 Transparent
- 558951 White(110/A1)
- 558954 White Opaque
- 558952 Yellow(210/A3)
- 558953 Brown(340/A4)



Catalyst – 1x2.5 gm

- 558955 High Viscosity Yellow(210/A3)
- 558956 Low Viscosity Yellow(210/A3)
- 558957 High Viscosity Transparent
- 558958 Low Viscosity Transparent
- 559224 Precision Applicator Tips

- 1 Total Etch®, 2 gm
- 1 Bottle Excite®, 5 gm
- 1 Monobond S®, 5 gm
- 1 Liquid Strip, 2.5 gm
- 5 Variolink®II Applicator Tips
- 1 Brush Holder
- 1 Vivapad®
- 1 Mixing Pad
- 2 VivaSticks®
- 50 Disposable Brushes
- Various Cannulas

Variolink II Try-In Standard Package - 558959

- 6 Syringes Variolink II Try-In; Shades Bleach XL(010), Transparent, White(110/A1), Yellow(210/A3), Brown(340/A4), White Opaque

Try-In – 1x2.5 gm

- 558960 Bleach XL(010)
- 548793 Transparent
- 548794 White(110/A1)
- 548795 White Opaque
- 548796 Yellow(210/A3)
- 548797 Brown(340/A4)

Intro Package - 560464

- 1 Syringe Variolink II Base 2.5 gm, Yellow(210/A3)
- 1 Syringe Variolink II Catalyst 2.5 gm; low viscosity; Transparent
- 1 Total Etch
- 5 Excite, 0.1 gm
- 1 Monobond S, 5 gm
- 1 Liquid Strip, 2.5 gm
- Various Accessories

Assortment - 558947

- 6 Syringes Variolink II Base 2.5 gm in shades Bleach XL(010), Transparent, White(110/A1), Yellow(210/A3), Brown(340/A4), White Opaque
- 2 Syringes Variolink II Catalyst 2.5 gm; high/low viscosity; Transparent
- 2 Syringes Variolink II Catalyst 2.5 gm; high/low viscosity; Yellow(210/A3)
- Various Accessories

Variolink Ultra Standard Package - 538319

- 1 Syringe Ultra Base, 3 gm, shade Yellow(210/A3)
- 1 Syringe Ultra Catalyst, 3 gm, shade Yellow(210/A3)

Multilink®

UNIVERSAL ADHESIVE CEMENT

Multilink System Pack

- 576862 Yellow Shade
- 576869 Transparent Shade
- 576864 White Opaque Shade



Multilink System Pack Contains:

- (1) Multilink Double Syringe, 6 gr
- (1) Multilink Primer A, 3 gr
- (1) Multilink Primer B, 3 gr
- (1) Multilink Metal primer, 5 gr
- Accessories - (50)brushes, (1)mixing pad, instructions for use

576825 Multilink Primer A, Primer B, contains 3 gr bottle

576827 Multilink Metal Primer, 5 gr

Multilink Refills - 1x6 gr

- 576865 Yellow Shade
- 576860 Transparent Shade
- 576866 White Opaque Shade



Vivaglass® CEM

GLASS IONOMER CEMENT

544414 Standard Pack

- 25 Capsules Vivaglass CEM

544413 Capsule Injector



USA

175 Pineview Drive, Amherst, NY 14228
716.691.0010 • 1.800.533.6825
716.691.2285 Fax
www.ivoclarvivadent.us.com

www.ivoclarvivadent.us.com

Call us toll free at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.
© 2004 Ivoclar Vivadent, Inc. Appeal, Multilink, Variolink and Vivaglass are trademarks of Ivoclar Vivadent, Inc.

CANADA

2785 Skymark Ave. Unit 1, Mississauga, Ontario L4W 4Y3
905.238.5700 • 1.800.263.8182 (Inside Canada)
905.238.5711 Fax
www.ivoclarvivadent.ca

100% CUSTOMER SATISFACTION
GUARANTEED!

ivoclar
vivadent
clinical